



Comfort
Release®

Global Biomedical
Technology

COMFORT RELEASE®

Transparent Dressings

INDICATIONS FOR USE

Use as a secondary dressing for partial-thickness wounds and full-thickness wounds.

Ideal for protecting post-surgical incisions, skin tears, lacerations, and abrasions.

Also helps prevent skin breakdown caused by friction.



ADHESIVE TECHNOLOGY

*Apply rubbing alcohol to the dressing border to deactivate the adhesive bond.



SKIN FRIENDLY

Designed for frail or sensitive skin, particularly for those with conditions or medications that compromise skin integrity.

TRANSPARENT WINDOW

- Waterproof
- Permits continuous observation
- Sterile barrier to external contaminants
- Moisture vapor permeable, allowing excess exudate to evaporate
- Flexes with the skin to promote comfort
- Low adhesive tack to secure tubing or primary dressing in place

ADHESIVE BORDER


- Low profile design to prevent edge rolling
- Secure adhesive that removes painlessly*



2.375" x 2.75"
GB 110 // HCPCS: A6257



5" x 5.25"
GB112-W // HCPCS: A6258



4.5" x 8.5"
GB 113 // HCPCS: A6258

Place Orders at: office@vmdinc.com // 561.588.3222

www.comfortrelease.com