



561-588-3222
office@vmdinc.com



USA

Manuka Honey

PRODUCT REFERENCE GUIDE

Naturally powerful healing
with the strength of Gentell.

gentell.com





We Make It Better.

At Gentell, we provide high-quality, affordable, patient-specific wound treatments to patients in nursing homes, hospices, wound care centers and other healthcare settings.

Founded in 1994 by President and CEO David Navazio, our history of innovation and growth has led to our position as the largest vertically-integrated wound care company in the world.

That same year, David worked with the National Association for the Support of Long-Term Care to advocate for having the Health Care Financing Administration recognize advanced wound care products. The result was a change to the Medicare Coverage Policies to recognize wound care and provide coverage for advanced wound care products for nursing home residents in the United States.

Our expansion into Canada, Brazil, Paraguay, India, the United Kingdom and New Zealand represents our ongoing commitment to delivering cost-effective and innovative wound care solutions to clinical customers and strategic partners around the world.

Today, Gentell is a global company with offices, distribution centers and manufacturing facilities on five continents. We set the standard for wound care products, education, customer service and advocacy. From our industry-leading advanced wound care products to our hands-on training, data-driven reporting and personalized consultations, we handle every aspect of wound care so you can focus on what you do best—caring for your patients.



MHCTGVMD_v1.1



CODE	PRODUCT	SIZE	PAGE
	About Manuka Honey		4
	What Makes Manuka Honey Unique?		5
	Why Honey for Wound Care?		5
MM0010	Manuka Honey Non-Adherent Dressing	4" x 5" (10cm x 12.5cm)	6
MM0011	Manuka Honey Non-Adherent Dressing	4" x 8" (10cm x 20cm)	6
MM0012	Manuka Honey Non-Adherent Wrap	4" x 48" (10cm x 122cm)	6
MM0013	Manuka Honey Fiber Premium	2" x 2" (5cm x 5cm)	7
MM0017	Manuka Honey Fiber Premium	4" x 5" (10cm x 12.5cm)	7
MM0030	Manuka Honey Fiber Select	2" x 2" (5cm x 5cm)	7
MM0031	Manuka Honey Fiber Select	4" x 5" (10cm x 12.5cm)	7
MM0070	Manuka Honey Fiber Dressing	2" x 2" (5cm x 5cm)	8
MM0071	Manuka Honey Fiber Dressing	4" x 5" (10cm x 12.5cm)	8
MM0048	Manuka Honey Wound Dressing	0.7oz (20g)	8
MM0049	Manuka Honey Wound Dressing	0.5oz (15g)	8
MM0051	Manuka Honey Wound Dressing	1.5oz (44g)	8
MM0053	Manuka Honey Skin Therapy Cream	1.75oz (50g)	9
	Manuka Honey Product Reference Sheet		10





Gentell's line of honey-based wound care dressings, gels and creams harness the proven power of medical-grade Manuka honey to promote healing, reduce inflammation and protect against infection.

Trusted by professionals worldwide, these products are backed by Gentell's industry-leading service, support and manufacturing capabilities, ensuring consistent quality and availability for your wound care needs.



Leptospermum scoparium (Manuka flower)

WOUND DRESSINGS

Gentell's honey-infused dressings create an optimal healing environment, managing moisture, protecting the wound bed and promoting faster recovery with the natural power of medical-grade Manuka honey.

GELS & CREAMS

Our honey-based gels and creams provide gentle, effective care for wounds requiring moisture balance and antimicrobial support. Ideal for sensitive skin, hard-to-dress areas and chronic wounds.




WHAT MAKES MANUKA HONEY UNIQUE?

Derived from the Manuka tree (*Leptospermum scoparium*), this type of honey is known for its unique antibacterial properties, attributed to compounds like methylglyoxal (MGO).

Our honey has the highest level of MGO of any medical grade honey available, with a minimum concentration of 560 mg/kg. This is important as suppurating wounds dilutes honey applied to wounds, reducing its effectiveness.

Manuka honey products offer consistent increased plasmin activity by inhibiting Plasminogen Activator Inhibitor (PAI), which digests the fibrin clot, loosening slough and necrotic tissue attached to the wound bed. Osmotic action draws lymph fluid into the wound, further assisting to loosen debris and necrotic tissue for easy removal.

All of our Manuka honey products are designed to support wound healing with natural ingredients with the strength of Gentell.



WHY HONEY FOR WOUND CARE?

Antibacterial Action - Honey contains compounds like hydrogen peroxide and has a high acidity and sugar content which inhibits bacterial growth and can help prevent wound infections.

Anti-Inflammatory Properties - Honey can reduce inflammation and swelling in the wound area to help manage pain and promote healing.

Antioxidant Activity - Honey is rich in antioxidants that help to neutralize harmful free radicals, protecting tissues and supporting the healing process.

Debridement - Honey can help remove dead tissue from the wound which is a crucial step in healing.

Moist Wound Healing - Honey can help maintain the moist wound environment that is important for optimal healing.



Honey has been used since ancient times and is understood by medical researchers to have beneficial uses in healing acute and chronic wounds.



Manuka Honey Non-Adherent Dressing

MM0010	4" x 5" (10cm x 12.5cm)	10/bx	20 bx/cs
MM0011	4" x 8" (10cm x 20cm)	10/bx	20 bx/cs

Manuka Honey Non-Adherent Dressing is an impregnated wound dressing with permeable acetate gauze with 100% active Leptospermum honey from New Zealand.

This dressing is indicated for use on skin grafts, surgical and arterial wounds, skin tears and burns with no to light (scant) exudate.

- High potency medical grade Manuka honey; pH of 3.0-4.5
- Reduces wound bed pH to optimize wound healing
- Promotes autolytic debridement and moist wound bed
- Manuka honey released into wound bed with gauze remaining on removal
- Flexible and comfortable



Manuka Honey Non-Adherent Wrap

MM0012	4" x 48" (10cm x 122cm)	4/box	6 bx/cs
---------------	----------------------------	-------	---------

Manuka Honey Acetate Gauze Wrap is an impregnated wound dressing wrap with permeable acetate gauze with 100% active Leptospermum honey from New Zealand.

This dressing is indicated for use on skin grafts, surgical and arterial wounds, skin tears and burns with no to light (scant) exudate.

- High potency medical grade Manuka honey; pH of 3.0-4.5
- Reduces wound bed pH to optimize wound healing
- Promotes autolytic debridement and moist wound bed
- Manuka honey released into wound bed with gauze remaining on removal
- Flexible and comfortable



Manuka Honey Fiber Premium

MM0013	2" x 2" (5cm x 5cm)	10/bx	20 bx/cs
MM0017	4" x 5" (10cm x 12.5cm)	10/bx	20 bx/cs

Manuka Honey Fiber Premium is an impregnated super absorbent gelling fiber pad for moderate to heavy exuding wounds. Flexible and comfortable, this dressing enables swelling and gelling when in contact with exudate and does not require a secondary absorptive dressing.

Indicated for moderate to heavy exudate on venous ulcers, diabetic foot ulcers, pressure ulcers, arterial and surgical wounds and skin tears.

- Wound dressing with super absorbent and cross-linked mesh fibers
- Highest medical grade Manuka honey concentration; most absorptive dressing with pH of 3.0-4.5
- Reduces or eliminates wound odor
- Reduces wound bed pH to optimize wound healing
- Stimulates autolytic debridement
- Eschar and debris tissue are taken up by the mesh fiber in the dressing resulting in a clean and viable wound bed



Manuka Honey Fiber Select

MM0030	2" x 2" (5cm x 5cm)	10/bx	20bx/cs
MM0031	4" x 5" (10cm x 12.5cm)	10/bx	20bx/cs

Manuka Honey Fiber Select is 100% *Leptospermum scoparium* honey from New Zealand impregnated into a pad of super absorbent gelling fiber for moderate to heavy exuding wounds.

Indicated for leg ulcers, pressure ulcers, first- and second-degree burns, diabetic foot ulcers, surgical wounds, trauma wounds and skin tears.

- Enables swelling and gelling when in contact with exudate
- Reduces wound bed pH to optimize wound healing
- Stimulates autolytic debridement
- Manuka honey is released into the wound bed with a sustained flow, while exudate and necrotic tissue are absorbed into the dressing fiber, resulting in a clean and viable wound bed
- Flexible and comfortable
- Does not require a secondary absorptive dressing



Manuka Honey Fiber Dressing

MM0070	2" x 2" (5cm x 5cm)	10/bx	20 bx/cs
MM0071	4" x 5" (10cm x 12.5cm)	10/bx	20 bx/cs

Manuka Honey Fiber Dressing is a wound dressing with super absorbent and cross-linked mesh fibers impregnated with Manuka honey and a pH of 3.0-4.5. Designed for ease of application and daily dressing changes for monitoring progress, these dressings can be used in combination with other products such as our Manuka Honey Wound Dressing for optimal medical grade honey dosage with every dressing change.

Indicated for moderate to heavy exudate on diabetic foot ulcers, pressure ulcers, venous leg ulcers and skin tears.

- Reduces wound bed pH to optimize wound healing
- Manuka honey released into wound bed with sustained flow, while exudate and necrotic tissue is absorbed into the dressing
- Eschar and debris tissue are taken up by the mesh fiber in the dressing resulting in a clean and viable wound bed
- Creates optimal wound bed healing environment and reduces or eliminates wound odor
- Reimbursable by insurance in the United States when wound meets coverage guidelines



Manuka Honey Wound Dressing

MM0048	0.7oz (20g)	1/bx	20bx/cs
MM0049	0.5oz (15g)	2/bx	50bx/cs
MM0051	1.5oz (44g)	1/bx	50bx/cs

Manuka Honey Wound Dressing contains 100% active Leptospermum honey from New Zealand in a tube applicator.

Indicated for use on venous ulcers, diabetic foot ulcers, pressure ulcers, arterial and surgical wounds, skin tears and second-degree burns with no to light (scant) exudate. May be used over-the-counter for minor abrasions, minor lacerations, minor cuts, minor scalds and burns.

- High potency medical grade Manuka honey with pH of 3.0-4.5
- Reduces wound bed pH to optimize wound healing
- Promotes autolytic debridement
- Manuka honey released into wound bed with sustained flow while exudate and necrotic tissue are absorbed into a dressing (applied separately)
- May assist with acne/hidradenitis suppurativa treatment



Manuka Honey Skin Therapy Cream

MM0053

1.75oz (50g)

1/ea

36/cs

Manuka Honey Skin Therapy Cream contains 30% Manuka honey for relief from chapped or cracked skin, eczema and burns.

This skin cream is the perfect solution for those who would benefit from a nourishing moisturizer and protection against burns in their day-to-day lives.

- Combines our medical-grade Manuka honey with soothing agents such as kukui oil and shea butter
- Restores skin moisture balance
- Protects and soothes dry, chapped skin
- Non-sticky, non-greasy, easy to apply
- Leaves a pleasant and long-lasting fragrance on the skin

ManukaMed is now Gentell Manuka Honey!

Use the product reference below to find your favorite product in its new packaging!

SKU	MANUKAMED NAME	GENTELL NAME	SIZE	QTY/BX
MM0010	MANUKAtex	Manuka Honey Non-Adherent Dressing	4" x 5" 10cm x 12.5cm	10
MM0011	MANUKAtex	Manuka Honey Non-Adherent Dressing	4" x 8" 10cm x 20cm	10
MM0012	MANUKAtex	Manuka Honey Non-Adherent Wrap	4" x 48" 10cm x 122cm	4
MM0013	MANUKAhd	Manuka Honey Fiber Premium	2" x 2" 5cm x 5cm	10
MM0017	MANUKAhd	Manuka Honey Fiber Premium	4" x 5" 10cm x 12.5cm	10
MM0030	MANUKAhd lite	Manuka Honey Fiber Select	2" x 2" 5cm x 5cm	10
MM0031	MANUKAhd lite	Manuka Honey Fiber Select	4" x 5" 10cm x 12.5cm	10
MM0048	MANUKApli	Manuka Honey Wound Dressing	0.7oz 20g	1
MM0049	MANUKApli	Manuka Honey Wound Dressing	0.5oz 15g	2
MM0051	MANUKApli	Manuka Honey Wound Dressing	1.5oz 44g	1
MM0053	MANUKAderm Skin Therapy Cream	Manuka Honey Skin Therapy Cream	1.75oz 50g	4
MM0070	MANUKAhd Superlite	Manuka Honey Fiber Dressing	2" x 2" 5cm x 5cm	10
MM0071	MANUKAhd Superlite	Manuka Honey Fiber Dressing	4" x 5" 10cm x 12.5cm	10
MM0121	MEDSAF	Fiber Dressing	4" x 5" 10cm x 12.5cm	10





Gentell is the world's largest vertically-integrated wound care company. From wound care, to nutritional services, to manufacturing and product design and development, Gentell's products and services make a positive impact on the lives of patients around the world.

VMD is proud to partner with Gentell to bring you these Manuka products. For more information or to place an order, please contact VMD at:

561.588.3222 | office@vmdinc.com