



Comfort
Release®

Global Biomedical
Technology

COMFORT RELEASE® Transparent Tape

INDICATIONS FOR USE

- Ideal for patients with sensitive or fragile skin, including the elderly, anyone on medications that can affect skin integrity, repetitive or chronic wounds, and post-op surgical patients.
- Recommended for those who frequently use adhesives or dislike the pain of tape removal.



FEATURES

- **Painless Removal** – Releases gently without causing trauma to the skin
- **Perforated** – Easy tearing in both directions—no scissors needed
- **Water-Resistant** – Stays secure even during bathing or showering
- **Breathable** – Perforations in the tape allow airflow to reduce moisture buildup and minimize maceration
- **Durable & Long-Lasting** – Stays in place to provide extended wear time
- **Activated Painless Release** – Designed for easy removal
- **Gentle on Skin** – Reduces risk of skin irritation or injury



ADHESIVE TECHNOLOGY

Apply rubbing alcohol to the tape's surface to deactivate the adhesive. Use an alcohol pad or soaked cotton ball, swiping the tape until it lifts from the skin.



WHY CHOOSE COMFORT RELEASE® TAPE?

- Cost-effective alternative to silicone tape
- Adheres stronger than silicone tape
- Transparent for easy skin monitoring
- Durable polypropylene for lasting performance

25 Rolls of Tape >
1" x 3.3' (2.5 cm x 100 cm)
GB107-05 // HCPCS: A4450



12 Rolls of Tape
1" x 4 Yards (2.5 cm x 365 cm)
GB103 // HCPCS: A4450



6 Rolls of Tape
2" x 4 Yards (5 cm x 365 cm)
GB104 // HCPCS: A4450

Place Orders at: office@vmdinc.com // 561.588.3222

www.comfortrelease.com