



DemeTECH[®]
"A Trusted Name in Medicine"

SURGICAL SUTURE

Portfolio

DEMETECH CORPORATION

DemeTECH Corporation
5980 Miami Lakes Dr,
Miami Lakes, FL 33014, USA.

Phone: (305) 824-1048

E-mail: info@demetech.us



www.demetech.us



MADE IN USA



Code Guide

Thread Type	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	
DemeSORB Polyglycolic Acid	PGA	DemeCRYL Polyglactin 910	G	DemeQUICK Rapid Absorbable	PGR	DemeGUT Chromic Catgut	CC	DemeGUTPLAIN Plain Catgut	CP	DemeCAPRONE Poliglecaprone 25	MO	DemeLON Nylon (Monofilament)	NL	DemeLON MULTI Nylon (Multifilament)	NM
DemeDIOX Polydioxanone	PX	DemeBOND Polyester (Braided)	PB	DemeDAC Polyester (Braided without covering)	PT	DemeLENE Polypropylene	PM	DemeFORCE Polyethylene(UHMPE)	FR	DemeSTEEL Surgical Steel	SW	DemeSILK Silk	SK	DemeDIOX BARBED Polydioxanone Barbed	PXB
DemeTECH PTFE Polytetrafluoroethylene	TE														

Thread Length	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	
45 cm	1	60 cm	6	70 cm	11	135 cm	16	25 cm	21	68 cm	26	35 cm	31	11 cm	36
75 cm	2	13 cm	7	120 cm	12	76 cm	17	300 cm	22	8 cm	27	6 cm	32	400 cm	37
100 cm	3	30 cm	8	180 cm	13	140 cm	18	500 cm	23	220 cm	28	4 cm	33	5.5 cm	38
125 cm	4	50 cm	9	20 cm	14	250 cm	19	40 cm	24	15 cm	29	3.5 cm	34	4.5 cm	39
150 cm	5	90 cm	10	200 cm	15	80 cm	20	185 cm	25	5 cm	30	10 cm	35	2.5 cm	40
														9 cm	41
														22 cm	42
														95 cm	43
														12 cm	44
														17.5 cm	45
														16 cm	46
														14 cm	47
														36 cm	48
														40 cm	49
														360 cm	50


Thread Color	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	
Silver	1	Beige	3	White	5	Colorless	7	Blue	9	Pink	11
Brown	2	Green	4	Black	6	Violet	8	Ivory	10	White/Blue	12

Thread Diameter	Code												
11/0	110	9/0	90	7/0	70	5/0	50	3/0	30	0	00	2	02
10/0	100	8/0	80	6/0	60	4/0	40	2/0	20	1	01	3	03
												4	04
												6	06

Needle Length	Code														
3,5	0,35	7	0,7	12	12	17	17	23	23	30	30	37	37	45	45
4	0,40	8	0,8	13	13	18	18	24	24	31	31	38	38	47	47
5	0,50	9	0,9	14	14	19	19	25	25	32	32	39	39	48	48
6	0,60	10	1,0	15	15	20	20	26	26	35	35	40	40	50	50
6,5	0,65	11	1,1	16	16	22	22	27	27	36	36	41	41	55	55
														58	58
														60	60
														65	65
														75	75
														76	76

Needle Curvature	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	
Straight (Double)	Ax2	1/2 Circle (Loop)	BL	1/4 Circle	D	1/4 Circle (Double)	E	3/8 Circle (Loop)	FL	5/8 Circle	H	5/8 Circle (Loop)	K	J Needle	J
Straight (Loop)	AL	1/2 Circle (Double)	B	1/4 Circle (Loop)	DL	3/8 Circle	F	3/8 Circle (Double)	G	5/8 Circle (Loop)	HL	Loop Round	I	Compound Curve (Double)	A
														Compound Curve	AL
														Ski	S

Cross Section	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	Code	
Round Bodied	0	Reverse Cutting	4	SBR Spatula	8	Reverse Cutting Prime	12	Curved Cutting SLIM BLADE	17	Curved Cutting Prime	21	Straight Cutting Cobra Point	25	Round Bodied Plus SLIM BLADE	29
Round Bodied Heavy	1	Reverse Cutting Heavy	5	Conventional Spatula	9	Precision Point Reverse Cutting	13	Round Bodied SLIM BLADE	18	CC	22	Round Bodied Plus	26	Round Bodied Super Heavy	30
Curved Cutting	2	Taper Cutting	6	CSU Spatula	10	Micropoint Reverse Cutting	14	Reverse Cutting SLIM BLADE	19	Blunt Point Heavy	23	RAFT Cutting	27	Deme Black VB	VB
Curved Cutting Heavy	3	Micropoint Curved Spatula	7	Precision Cosmetic Standard Prime	11	Blunt Point	16	Taper Cutting Heavy	20	OverStitch	24	Round Bodied Plus Heavy	28		



Needle Chart

Round Bodied

1/2 Circle

DRB-2		13mm
DRB-1		17mm
DSH-1		22mm
DSH		26mm

1/2 Circle (LOOP)

DCT LOOP		40mm
----------	---	------

3/8 Circle

DC-1		13mm
DBB		17mm
DTE		31mm

3/8 Circle (DOUBLE)

DC-1		13mm
------	---	------

5/8 Circle

DUR		26mm
-----	---	------

Taper Cutting

1/2 Circle

DV-5		17mm
DV-6		22mm
DV-34		36mm
DV-40		48mm

1/2 Circle (LOOP)

DV-7		26mm
------	---	------

3/8 Circle

DV-4		17mm
------	---	------

J Needle

DV-38		31mm
-------	---	------

Precision Point Reverse Cutting

1/2 Circle

DPS-4		16mm
-------	---	------

1/2 Circle

DPS-3		16mm
DPS-2		19mm
DPS-1		24mm
DPS		26mm
DPSL		30mm

Reverse Cutting

1/2 Circle

DJ-1		17mm
DC-1		22mm
DCP-2		26mm
DCP-1		36mm

3/8 Circle

DC-3		13mm
DFS-3		16mm
DFS-2		19mm
DFS-1		24mm

3/8 Circle (DOUBLE)

DLR		75mm
-----	---	------

Straight

DKS		60mm
-----	---	------

Curved Cutting

3/8 Circle

DCFS-3		16mm
DSB-25		26mm
DCFSLX		40mm

Straight

DKS		60mm
-----	---	------

Micropoint Curved Spatula

3/8 Circle

DV-4		17mm
------	---	------

3/8 Circle (DOUBLE)

DV-38		31mm
-------	---	------

Precision Cosmetic Standard Prime

3/8 Circle

DPC-1		13mm
DPC-3		16mm
DPC-5		19mm
DPC-26		26mm

| Surgical Sutures

DemeTECH Corporation is a world-renowned leader in surgical sutures, blades, and staplers. All DemeTECH suture lines are manufactured with the highest quality medical grade materials.

Absorbables Sutures

DemeCAPRONE

Poliglecaprone 25 Suture

DemeSORB

Polyglycolic Acid Suture

DemeQUICK

Rapid Absorbable Suture

DemeCRYL

Polyglactin 910 Suture

DemeGUT

Chronic Catgut Suture

DemeGUT PLAIN

Plain Catgut Suture

DemeDIOX

Polydioxanone Suture



Non-Absorbable Sutures

DemeFORCE

Polyethylene(UHMWPE)

DemeBOND

Polyester Suture(Braided)

DemeLON

Nylon Suture(Monofilament)

DemeLON MULTI

Nylon Suture (Multifilament)

DemeLENE

Polypropylene Suture

DemeDAC

Polyester Suture (Braided without covering)

DemeSILK

Silk Suture

DemeSTEEL

Surgical Steel Suture


DemeDIOX BARBED

Polydioxanone Barbed Suture

DemeTECH PTFE

Polytetrafluoroethylene

| Primary Label



- 1- Metric
- 2- Thread Length
- 3- UDI Barcode
- 4- Lot
- 5- Expiration Date
- 6- Sterilization Method
- 7- USP
- 8- Needle Cross Section
- 9- Needle Curvature
- 10- Sutures Material
- 11- Thread Description
- 12- Needle Shape
- 13- Suture Code