




167 West Orangethrope Avenue, Placentia, CA 92870



SCAN ME

Duraholder® IPS



Biosealnet.com

To Learn More
Please Contact:




ISO 13485

Our unique fabrication capabilities allow us to offer a wide range of applications. Please note Duraholder® IPS is not intended to provide a sterile barrier and should be used in conjunction with an outer wrap or container system.


Note: Please allow for a 1” tolerance on outer dimensions around the sealed pockets. More duraholders are available in various sizes please see website for more.

Part Number	Description	Size	UOM
A DURA050/100	Duraholder 1 Pocket	10" x 9 ¾"	100/cs
B DURA100/100	Duraholder 2 Rows, 6 Pockets	11 ½" x 17 ½"	100/cs
C DURA102/100	Duraholder 2 Rows, 6 Pockets	11 ½" x 17 ½"	100/cs
D DURA200/100	Duraholder 2 Rows, 6 Pockets w/ Notches, w/ Tape	10 ½" x 18 ½"	100/cs


A Duraholder 50




B Duraholder 100



C Duraholder 102




D Duraholder 200




Part Number	Description	Size	UOM
E DURA300/100	Duraholder 2 Rows, 6 Pockets	18 ½" x 16"	100/cs
F DURA317/100	Duraholder 1 Row, 18 Pockets	19 ½" x 11"	100/cs
G DURA415/100	Duraholder 1 Pocket	4 ⅝" x 10 ⅝"	100/cs
H DURA505/100	Duraholder 1 Pocket	5 ½" x 5 ¼"	100/cs
I DURA512/100	Duraholder 1 Pocket	12" x 5"	100/cs
J DURA602/100	Duraholder 2 Rows, 6 Pockets w/ Slits	18" x 16"	100/cs


E Duraholder 300




F Duraholder 317




G Duraholder 415




H Duraholder 505



I Duraholder 512



J Duraholder 602



The Duraholder® Instrument Protection System (IPS) was designed to prevent costly damage to delicate neuro instruments and eliminate peel packs. Duraholder® is made using HALYARD® ONE-STEP Sterilization Wrap, adhering to AAMI guidelines.

Protect Every year, hospitals spend millions of dollars on instrument repair and replacement; money which could be better spent serving patients. Bioseal's Duraholder® is a proven method of protecting delicate instruments from the sterilization facility all the way to the operating room floor – reducing instrument repair costs.

Organize On average, it costs up to \$65 for every minute spent in the O.R. Disorganized instrument sets can add to the already expensive time constraints, costing hospitals thousands of dollars. Duraholder®, an inventive Instrument Protection System designed specifically to help you keep track of the instruments, resulting in fewer lost instruments, and reduced time spent organizing instruments in the O.R.

Comply Duraholder® adheres to sterility and sterile product guidelines and protocols set by AAMI and ASTM. The AAMI guideline ST79 states that “paper-plastic pouches are not appropriate for use within wrapped sets or containment devices”. The Duraholder® fits nicely into the new AAMI standards.

Take a look at our most popular configurations and see why our Duraholder® products have become a customer favorite!