

GOOD WORKS®

PERFORMANCE PRO®
NITRILE EXAM GLOVES



BOSMA Enterprises®

Good Works® Performance Pro® examination gloves from Bosma Enterprises offer textured fingers for a strong grip and enhanced tactile sensitivity. The comfortable fit and feel makes this glove great for use in a variety of settings. Every purchase from Bosma creates jobs for people who are blind or visually impaired.

BOA#: 36C10X23G0002

Packaging: 100/Box; 10 Boxes/Case

Color: Blue

Finish: Textured fingertips

Type: Ambidextrous/Non-Sterile /Disposable/Single Use

Length: 9.45 inches (min)

Cuff: Beaded

Palm Thickness: 3.5 mil

Finger Thickness: 4.5 mil

Cuff Thickness: 3 mil

Tensile Strength: 14 Mpa (min)

Elasticity/Elongation: 500% (min)

Quality Standards:

- Tested for use with chemotherapy drugs in accordance with ASTM D6978
- Complies with FDA Biocompatibility Guidance for Medical Devices
- Manufactured to ASTM D6319 standard for nitrile gloves for medical application
- ISO 9001 and 13485
- 510(K) / Approved for medical use
- Component materials comply with FDA regulations for food contact
- Manufactured by American Nitrile in Grove City, OH.



Bosma #	NSN	Size
72-50601	6515-00-NIB-9762	Small
72-50602	6515-00-NIB-9763	Medium
72-50603	6515-00-NIB-9764	Large
72-50604	6515-00-NIB-9765	X-Large

For more information go to bosma.org or call **800.362.5463**

BOSMA Enterprises®

6270 Corporate Drive
Indianapolis, IN 46278

Bosma Enterprises®, creating opportunities for individuals who are blind or visually impaired. Your purchase creates employment for more than 100 people who are blind.



ISO 9001:2015
ISO 13485:2016
CERTIFIED