

# GRANITE™

POWDER-FREE NITRILE DISPOSABLE GLOVES



AMERICAN NITRILE

STRENGTH AND SENSITIVITY AT THE TIP OF YOUR FINGERS



- ✓ Chemical Resistant Nitrile
- ✓ Optimal Comfort and Protection
- ✓ Enhanced Dexterity
- ✓ Greater Tactile Sensitivity
- ✓ Latex Free

Learn more or inquire at:  
[americannitrile.com](http://americannitrile.com)

# GRANITE™

## POWDER-FREE NITRILE DISPOSABLE GLOVES



### AMERICAN NITRILE

#### AMERICAN NITRILE ADVANTAGES

Every American Nitrile glove is proudly manufactured at our state-of-the-art facility in Grove City, Ohio. As the preferred domestic source for nitrile gloves, we ensure:

- ✓ **Delivery:** Stress-free and just-in-time, eliminating unpredictable trans-pacific shipping, port delays, and the need to carry high safety stock.
- ✓ **Corporate social responsibility:** Sustainability in energy and raw material usage plus visibility into ethical labor practices.
- ✓ **Product quality:** Our gloves are 100% made in the USA to ensure safety and quality are never compromised.

#### KEY FEATURES AND BENEFITS

##### Tested

Meets US FDA device license regulations, complies with ASTM standards, and cleared for use with chemotherapy drugs in accordance with ASTM D6978.

##### Authentic

US-based production and delivery ensures quality and instills confidence in the safety and cleanliness of our products.

##### Wearability

Smooth inner surface for easy donning.

##### Superior AQL

Our gloves meet a pinhole AQL of 1.5, 60% better than the FDA requirement of 2.5.

#### GLOVE SIZE

#### REORDER NUMBER

Small	ANGR552
Medium	ANGR553
Large	ANGR554
X-Large	ANGR555
XX-Large	ANGR556

#### PHYSICAL PROPERTIES

	Typical Performance		Testing Method
Length	245 mm/9"		ASTM D3767
Freedom from Holes	2.5 AQL		ASTM D5151
Palm Thickness	0.135 mm/5.3 mil		ASTM D3767
Finger Thickness	0.17 mm/6.6 mil		ASTM D3767
Glove Weight Medium	5.5 grams		
	Before Aging	After Aging	
Ultimate Tensile Strength (MPa)	>27	>24	ASTM D412 & D573
Elongation at Break (%)	>500	>400	ASTM D412

#### PRODUCT INFORMATION

Material	Nitrile
Color	Black
Glove Design	Chlorinated, Powder Free
Cuff	Beaded
Audit Standards	GMP
Packaging	100 Gloves Per Box, 1000 Gloves Per Case
Shelf Life	3 years
Country of Origin	USA
User Needs	Chemical Splash and Abrasion Protection, Puncture Resistance
Sizes	S, M, L, XL, XXL

#### SUSTAINABILITY IMPACT:



Low carbon footprint, no transpacific vessels used



Up to 40 gallons of water recycled per case



100% recycled/recyclable cardboard packaging



Powder-Free



Fingertip Texture



Ambidextrous



Non-Sterile



Single Use Only



510(k)  
Cleared  
Products