

CLEANPatch[®]-V

PRODUCT INFORMATION*

Item Number	CleanPatch [®] -V Size	Application
#21001 (10/box)	1.75 x 4 inch strip (S)	Damage up to 2 inches
#22001 (5/box)	5.75 x 8 inch rectangle (R)	Damage up to 6 inches
#23001 (10/box)	7 S + 3 R	Damage up to 6 inches
#24001 (3/box)	2 S + 1 R	Damage up to 6 inches

* SKU availability may vary by region. Please contact your local distributor for more information.



Sanitis
Austin, TX
(888) 789-6266
www.SanitisGlobal.com

Manufactured by
Surface Medical Inc.
Calgary, AB, Canada T2S 1W2
info@surfacemedical.ca

© Surface Medical Inc. 2019

CPV09182019

CLEANPatch[®]-V

*Medical Surface Repair patch for vinyl upholstery.
Restores damaged surfaces to an intact and hygienic state.*



WHY CLEANPATCH[®]-V?

Vinyl equipment is easily damaged from daily wear and tear, and it is important to repair damage quickly. When vinyl is damaged, it cannot be properly disinfected and may harbor dangerous pathogens, increasing the risk of HAIs among patients.

Until now, there was no safe solution for repairing damaged vinyl. CleanPatch[®]-V is designed to repair damage quickly, easily and safely — and for a fraction of the cost of replacement.

BEFORE CLEANPATCH®-V

Replace or Recover Equipment	✓ Validated ✗ Unaffordable ✗ Inconvenient
Tape Over Damage	✓ Inexpensive ✗ Not Validated ✗ Not Durable
Do Nothing	✗ Risk Patient Safety




THE CLEANPATCH®-V SOLUTION

CleanPatch®-V is the first medical surface repair product for vinyl equipment such as wheelchairs, examination tables and geriatric chairs. It restores damaged surfaces to an intact and hygienic state — ready for cleaning!

Clinically Tested	<ul style="list-style-type: none"> Fully cleanable with hospital disinfectants Adopted by healthcare facilities globally FDA, Health Canada, CE Mark, and TGA registered
Safe and Simple	<ul style="list-style-type: none"> Early intervention reduces the risk of cross contamination Impervious to fluids under rigorous conditions Simple peel and stick application, adheres in seconds
Cost Effective	<ul style="list-style-type: none"> Extends the life of vinyl equipment Improves patient safety Significant savings compared to equipment replacement

CleanPatch®-V should be used as part of a regular equipment maintenance program

	✓ Use CleanPatch-V	✗ Do Not Use CleanPatch-V
Equipment Condition	Within service life of equipment No signs of fluid ingress	No longer suitable for use Surface is stained or warped
Damage Size and Location	Less than 6 inches (15 cm) Flat areas or over edges (2-D) Not over seams	More than 6 inches (15 cm) Sharp corners (3-D) Over seam
Application	One patch per damaged area	Layering multiple patches is not recommended
Examples	 	 

To ensure product integrity, do not cut CleanPatch®-V



References:

FDA Safety Communication, April 2013
 Oliveira, AC. et al. Journal of Nursing UFPE on line, Jan 2013
 Otter JA, et al. Infection Control and Epidemiology, Vol 32 (7) July 2011
 Viana, REI H. et al. American Journal of Infection Control, 1;44(4) April 2015